



# PHAROS™ HD

Revolution in Headlight Imaging



# PHAROS™ HD

## Revolution in Headlight Imaging

Pharos™ HD is the most advanced coaxial headlight imaging system today. Providing video and still image capability in full 1080p/60 HD, it combines professional-grade single-chip CMOS HD technology integrated with our Hatteras™ LED fiber optic light source system. Compatible with most external HD recording devices, Pharos™ HD is lightweight, comfortable, easy to use and unequalled for surgeon's point-of-view image capture.

\* Available with our Classic or Key West™ "Bikini" headbands.

### Pharos™ Specifications

- 1080p/60 Video Output
- Camera Chip – CMOS Digital
- 1920 x 1080 True HD
- Image Stabilization
- 1 SDI output, 1 HDMI output
- White Balance, Brightness and Zoom Functions
- Push-Button Membrane Controls
- 4.5" Diameter spot at 16" working distance (11.4cm @ 40cm)
- 1-Year Limited Manufacturer's Warranty

### System Includes

- Pharos™ HD Headlight Camera Unit with 10' (Approx. 3m) Fiber Optic Cable
- Hatteras™ LED High-Intensity Light Source
- Pharos™ HD Camera Control Unit
- Appropriate Video Cables and Power Cords
- Mobile Stand

### Product Codes

- BFW 8010 SYSTEM – Pharos™ HD Headlight Imaging System – Classic Headband
- BFW 8012 SYSTEM – Pharos™ HD headlight Imaging System – Key West "Bikini" Headband



Powered by:

Long Island  
**Technology**  
Group