

EXAMVISION™

KEPLER
ADVANCED
By ExamVision

FEEL
MASTERY
WEAR
POWER



Get the perfect magnification level for all your procedures with only one loupe.



Turn N'Click system

Kepler Advanced is the groundbreaking adjustable magnification loupe brought to you by ExamVision. Thanks to its *Turn N'Click* system, switch magnifications easily in a single motion to achieve perfect precision for any procedure.

Enjoy the following benefits that make Kepler Advanced rise above the adjustable magnification loupe market:



4 magnifications in 1 loupe

The only loupe with 4 magnifications, for unrivalled versatility. Magnifications: 3.6x-4.5x-5.5x-6.4x



Highest adjustable magnification on the market

Up to 6.4x magnification, for enhanced precision and ergonomics



The world's lightest adjustable loupe

Kepler Advanced is more than 20 grams lighter than other comparable loupes, for ultimate comfort



Durable quality components

Handmade in Denmark with titanium and carbon fibre composite, and featuring world-class, cutting-edge lenses

The world's lightest adjustable magnification loupe, fully loaded with exceptional features



1. Four magnifications in one loupe, with easy adjustment via its *Turn'N'Click* system: 3.6x - 4.5x - 5.5x - 6.4x

2. Lightweight from just 75 grams

3. Adjustable declination angle up to 40° for optimised posture and tailored working distance from 32 cm to 65 cm

4. Unique advanced prismatic system

5. Made with titanium and carbon fibre reinforced composite

6. Optional BlueLight Protection



7. Compatible with ExamVision's LED light systems

8. Optional Eye Protect System™ available

9. Accommodates your prescription



Advanced optical technology for all your magnification needs

From diagnostics to near-microscopic magnification needs, Kepler Advanced ensures you achieve perfect sharpness of vision, while **maintaining a comfortable, ergonomic position.**

Kepler Advanced is **specifically designed for expert users within the dental and medical sectors who want to upgrade to a more powerful, high-end loupe.** It is intricately handmade with exclusive loupe technology and sophisticated materials.

At ExamVision we explore and master **optical technologies to offer you superior solutions.** Three Red Dot Design Awards are testament to this commitment to excellence. Our loupes are 100% custom-made for all your work and vision demands.

Get in touch with our highly trained specialists and see for yourself how far Kepler Advanced will take you.


→ Key magnification specifications:

Adjustable magnifications <i>Easy Turn 'N' Click Technology</i>	3.6x	4.5x	5.5x	6.4x	Working distance
Field of view (at 40 cm)	7 cm	6 cm	5 cm	4 cm	From 32 to 65 cm
Depth of Field (at age 40 years)*	8 cm	6 cm	5 cm	3.5 cm	

* Can vary with accommodation ability. Younger eyes will achieve greater depth of field.

→ Kepler Advanced is ideal for the following:

- All-day use for general procedures and diagnoses
- Intricate procedures that require high magnification, while providing the right depth of field for optimised ergonomics
- Near-microscopic magnification with the freedom and comfort of a non-static position





KEPLER
ADVANCED
By ExamVision

**Your upgrade awaits.
Ready for mastery and power?**
Contact your local dealer for your very
own pair of Kepler Advanced loupes.

www.exam-vision.com

Follow us on:

 facebook.com/ExamVisionDK

 instagram.com/examvision/