


# PUREVUE™ VISUALIZATION SYSTEM

Differentiation That Matters, So You Can **Visualize Success**



# PUREVUE™ Visualization System

A compact, scalable surgical imaging system delivering excellent image quality, brilliant live video and comprehensive image management.



**VISUALIZE  
SUCCESS**



## Excellent Image Quality

- Proprietary, purpose-built endoscopy camera sensor delivers real-time, native 1080P HD video images without upscaling, for a crystal clear view of what's happening, when it's happening
- Autoclavable camera featuring full HD CMOS technology designed exclusively for endoscopy
- Ergonomic and intuitive design with responsive tactile buttons
- High quality display options, from Value HD to 4K
- In-joint visibility enhanced when used in combination with our VAPR VUE® Radiofrequency System and FMS VUE™ Fluid Management and Tissue Debridement System



## Image & Case Management

- Comprehensive, easy-to-set-up Image and Case Management
- Access and control remotely via secure mobile app
- Aids post-op surgeon-to-patient communication
- Integration to facility record archives



## Compact, Easy-To-Use & Scalable

- Ease of use and reduced footprint with the camera control, light source and local image management merged into one unit
- Intuitive user interface with streamlined menu structure
- Fast picture capture display, for confidence you got the shot you want without delay
- Built-in forward compatibility with future software and hardware
- Flexible configurations to meet facility needs now and into the future

## Worry Free Ownership

- Full range of optical accessories to meet your needs
- Comprehensive service and repair offerings
- Part of the **DePuy Synthes VUE™ Arthroscopic Equipment Solutions** portfolio



VAPR VUE® Radiofrequency System

FMS VUE™ Fluid Management and  
Tissue Debridement System





## EXCELLENT IMAGE QUALITY

Crystal clear, real-time video



## COMPACT, EASY-TO-USE & SCALABLE

Small footprint, streamlined interface, forward compatibility



## COMPREHENSIVE IMAGE & CASE MANAGEMENT

Robust, yet simple and intuitive



## WORRY FREE OWNERSHIP

Delivering accessory breadth and peace of mind

## Visualization Components & Accessories

CODE	DESCRIPTION
242302	Camera Control System with Light Source (ROW)
242401	Camera Head with C-Mount, Autoclavable
242283	Video Cart with 230V Integrated Power Source
242284	Primary Multi-directional Display Mount
242518	Secondary Multi-directional Display Mount
242519	Rotatable Tablet Arm
242285	Secondary Display Cart
242436	Coupler, Zoom - Autoclavable
242311	Coupler, C-Mount, 24mm, Cold Sterilized
242312	Coupler, C-Mount, 17mm, Cold Sterilized
242316	Coupler, C-Mount 20mm, Cold Sterilized
242671	Medical Display, NDS Radiance Ultra 32" HD
242672	Medical Display, NDS Radiance Ultra 32" 4K/UHD
242351	Medical Display, NDS OHD2 26" High Bright Surgical
242325	Image Management: Photo and Video Capture System with Wireless Adaptor

CODE	DESCRIPTION
242370	Sony UP-DR80MD Printer
242349	Sony A4 Ink & Paper RF UPC-R80MD
242371	Sony 8.5 x 11 Ink & Paper RF UPC-R81MD
242386	Light Guide, OLYM-ACMI, 3.5MMX3.0M
242389	Light Guide, UNIV, 3.5MMX3.0M
242391	Light Guide, OLYM-ACMI, 3.5MMX4.0M
242390	Light Guide, UNIV, 5.0MMX3.0M
242387	Light Guide, OLYM-ACMI, 5.0MMX3.0M
242388	Light Guide, OLYM-ACMI, 5.0MMX4.0M
242335	Visualization Tray (Square)
242336	Visualization Tray Lid (Square)
242337	Visualization Tray (Scope, Rectangular)
242338	Visualization Tray Lid (Scope, Rectangular)
242356	ZeroWire® Duo, Tx/Rx Set, DVI + SDI Input, 3-Channel
242357	ZeroWire® Accessory Kit, for 24"/26" display

\* Additional video cart power cords, video displays and display cables available. Speak to your DePuy Synthes sales consultant for more information.

## Scopes, Sheaths, Obturators and Trocars

SCOPE TYPE	SCOPE SIZE	SHEATH		OBTURATOR	TROCAR
		SINGLE STOPCOCK	DUAL STOPCOCK		
STORZ LOCK	(242047) 1.9, 0 deg, 60 mm	242561	242523	242556	242555
	(242048) 1.9, 30 deg, 60 mm	242562	242524	242556	242555
	(242055) 2.7, 0 deg, 75 mm	242565	242528	242560	242559
	(242056) 2.7, 30 deg, 75 mm	242566	242529	242560	242559
	(242058) 2.7, 70 deg, 75 mm	242567	242540	242560	242559
	(242005) 4.0, 0 deg, 175 mm	-	242547	242585	242584
	(242006) 4.0, 30 deg, 175 mm	-	242548	242585	242584
	(242008) 4.0, 45 deg, 175 mm	-	242550	242585	242584
	(242009) 4.0, 70 deg, 175 mm	-	242551	242585	242584



PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

**DePuy Mitek, Inc.**  
325 Paramount Drive  
Raynham, MA 02767  
T. +1 (800) 382-4682

**Medos International SÁRL**  
Chemin-Blanc 38  
2400 Le Locie  
Switzerland

Distributed in the USA by:  
**DePuy Mitek, Inc.**  
325 Paramount Drive  
Raynham, MA 02767  
USA  
Tel: +1 (800) 382-4682

Authorized European Representative:  
**DePuy International, Ltd.**  
St Anthony's Road  
Leeds LS11 8DT  
England  
Tel: +44 (0) 113 387 7800



[www.depuyorthos.com](http://www.depuyorthos.com)