

# ULTRASEAL™


FEMORAL PRESSURISERS



# Pressurisation - The critical step in modern cementing technique<sup>1</sup>

Interdigitation of cement into cancellous bone is essential for long term fixation of the femoral stem. To achieve it requires optimal cementing technique using the best cementing accessories. DePuy Femoral Pressurisers have demonstrated their effectiveness in this demanding area.


## Femoral Cement Pressurisation



### Standard Femoral Pressuriser

Flexible silicon base seals the femur.

Backplate allows maximum pressurisation during cement insertion.



### Wedge Femoral Pressuriser

Strong simple shape fits and seals the femur (with the greater trochanter removed) allowing pressurisation during cement insertion.

1. Herberts, Malchau, Gothenburg, Prognosis of Total Hip Replacement 1996.
2. Independent testing - Unit for Joint Reconstruction - Oswestry. UK 1999.

#### Ordering Information

3206-002	Wedge Femoral Pressuriser	(pack of 5)
3206-005	Standard Femoral Pressuriser	(pack of 5)

This publication is not intended for distribution in the USA

CMW® is a registered trademark and Ultraseal™ is a trademark of DePuy International Ltd.  
© 2004 DePuy International Ltd. All rights reserved.

Cat No: 4010-070



**DePuy International Ltd**  
Trading as: **DePuy CMW**  
Cornford Road  
Blackpool FY4 4QQ  
England  
Tel: +44 (1253) 765 167  
Fax: +44 (1253) 697 431

